


(a) Procedure

1. Using a 50 cm³ measuring cylinder, transfer 30 cm³ of I₂(aq) into a 250 cm³ beaker. To the same beaker, add 30 cm³ of H₂SO₄(aq) using another 50 cm³ measuring cylinder.
2. Measure 30 cm³ of CH₃COCH₃(aq) into a 50 cm³ measuring cylinder.
3. Add the CH₃COCH₃(aq) into the same 250 cm³ beaker in step 1. Start the stopwatch immediately. Stir the mixture using a glass rod to ensure even mixing.
4. Using a 10.0 cm³ pipette, transfer 10.0 cm³ of the reaction mixture into a 250 cm³ conical flask. At time 2 min, add 10 cm³ of NaHCO₃(aq), measured using a 10 cm³ measuring cylinder, into the 250 cm³ conical flask.
5. Immediately titrate the iodine present in the conical flask against Na₂S₂O₃(aq) from a burette, until the solution turns pale yellow. Add 1 cm³ starch from a dropper and continue the titration until the blue-black solution in the conical flask turns colourless.
6. Repeat steps 4 and 5 at 6, 10, 14, 18, 22, 26 minutes.
7. Record the titration results in a suitable table.

(b)



The volume of Na₂S₂O₃ used for titrating a sample of quenched aliquot at time *t* is proportional to the concentration of I₂ present in the reaction mixture at time *t*.

Since the gradient of the graph is constant with a negative value, the rate of reaction is constant. Hence, the order of reaction with respect to iodine is zero since the reaction rate is independent of the concentration of iodine.