Anglo-Chinese Junior College

JC2 Biology Preliminary Examination Higher 2

CANDIDATE NAME		FORM CLASS	
TUTORIAL CLASS		INDEX NUMBER	

BIOLOGY

Paper 4 Practical

Candidates answer on the Question Paper. Additional Materials: As listed in the Confidential Instructions.

READ THESE INSTRUCTIONS FIRST

Write your Name, Class and Index number in the spaces on all the work you hand in. Give details of the practical shift and laboratory, where appropriate in the boxes provided. Write in dark blue or black pen.

You may use an HB pencil for any diagrams or graphs.

Do not use staples, paper clips, glue or correction fluid.

Answer all questions in the spaces provided on the Question Paper.

The use of an approved scientific calculator is expected, where appropriate. You may lose marks if you do not show your working or if you do not use appropriate units.

At the end of the examination, fasten all your work securely together. The number of marks is given in brackets [] at the end of each question or part question.

For Examiners' use		
1	/ 12	
2	/ 23	
3	/ 20	
Total	/ 55	







9744/04

06 August 2024 2 hours 30 minutes

Answer **all** questions.

1 In some countries, certain plants have seasonal growth. During the summer plants transport sucrose from the leaves to store it as starch in the roots.

Table 1.1 summarises the changes in the leaves and roots during the year.

Table 1.1

season	leaves	roots
summer	leaves synthesise sucrose	sucrose stored as starch
winter	no leaves	starch stored
spring	leaves growing	starch converted to glucose

You are required to identify the source of four plant extracts. These have been taken from

- a root in winter
- a root in spring
- phloem sap from the stem in summer
- phloem sap from the stem in winter
- (a) Use the information in Table 1.1 to predict which substances you would expect to be present in each of the four plant extracts, then complete Table 1.2.

Key: \checkmark (tick) substance present in plant extract

X (cross) substance absent from plant extract

Table	1.2
-------	-----

source of plant extract	substances present in each of the plant extracts			
	starch	sucrose	glucose	
root in winter	\checkmark	X	X	
root in spring	√ or X	X	\checkmark	
phloem sap from the stem in summer	X	\checkmark	X	
phloem sap from the stem in winter	X	X	X	

[2]

1m for first two rows, 1m for last two rows;

One sample was taken from each of the four plant extracts shown in Table 1.2.

You are required to identify from which plant extract each of the four samples **S1**, **S2**, **S3** and **S4** was taken.

You are provided with:

- Benedict's solution, in a container labelled B
- iodine solution, in a container labelled I
- hydrochloric acid, in a container labelled H
- sodium hydrogencarbonate powder, in a container labelled S
- a source of heated water.

Hydrochloric acid H and sodium hydrogencarbonate S are irritants, Benedict's solution B is harmful and iodine solution I is a stain. Suitable eye protection should be worn. If any of these reagents come into contact with your skin, wash off immediately under cold water.

Do not carry out any tests until you have read the instructions on pages 3 to 5.

(b) Describe the two tests that show that starch and glucose are present in a plant extract.

test for starch:

On a <u>white tile,</u> place <mark>3 drops</mark> of <u>sample</u>, then add <mark>1 drop</mark> of <u>iodine;</u> OR To <mark>1 cm³</mark> of <u>sample,</u> add <mark>3 drop</mark> of <u>iodine;</u>

.....[1]

test for glucose:

To <u>2 cm³</u> of <u>sample</u>, add <u>2 cm³</u> of <u>Benedict's solution</u> (volume must be 1:1);
 Heat in a <u>boiling water bath</u> for <u>2 minutes</u> (min 1 minute, max 3 minutes);

......[2]

For the test for sucrose, refer to steps 1 to 9.

1 Before starting the test for sucrose, perform the test for glucose you have described in (b).

4

- (c) Explain why the test for glucose must be conducted before the test for sucrose can be carried out.
 - 1. The test for sucrose / non-reducing sugar test involves <u>hydrolysis</u> of sucrose / non-reducing sugar to form glucose / reducing sugar;
 - Hence necessary to test for the absence of glucose / reducing sugar at the start; OR
 - 3. The test for sucrose / non-reducing sugar will eventually still involve a test for glucose / reducing sugar / Benedict's test;
 - 4. Hence necessary to test for the absence of glucose / reducing sugar at the start; [2]
- **2** To start the test for sucrose, put 2.0 cm³ of a fresh sample into a test-tube.
- **3** Put 2.0 cm³ of **H** into the same test-tube.
- 4 Shake the test-tube gently to mix the contents.
- **5** Put the test-tube in a boiling water-bath. Leave the test-tube for 2 minutes.
- 6 After 2 minutes, remove the test-tube from the water-bath and put it in a test-tube rack.
- 7 Leave the test-tube to cool for a further 3 minutes.
- 8 After 3 minutes, put a small amount of **S** into the test-tube. The mixture will fizz and rise up the test-tube. Continue to add small amounts of **S** until there is no more fizzing.

Note: there may be a little of S left in the bottom of the test-tube. This will not affect the results.

9 Perform the test for glucose you have described in (b) on the solution from step 8.

Proceed as follows.

Use the beaker labelled **hot water** to collect approximately 400 cm³ of hot water from where it is provided in the laboratory. Heat the water to boiling, if needed.

Select the appropriate reagents from those provided and carry out suitable tests to identify the samples **S1**, **S2**, **S3** and **S4**.

(d) Record your observations in an appropriate table. You do **not** need to make conclusions on the presence or absence of each substance in the samples tested.

	Observations			
Sample	Test for starch / lodine test	Test for glucose / Reducing sugar test	Test for sucrose / Non-reducing sugar test	
S1	The mixture remained yellow-brown	The mixture remained blue	The mixture remained blue	
S2	The mixture turned blue-black	The mixture remained blue	The mixture remained blue	
S3	The mixture remained yellow-brown	The mixture remained blue	Brick-red ppt is formed	
S4	The mixture turned blue-black	Brick-red ppt is formed	Brown/Brick-red ppt is formed	

- 1. Appropriate table layout (sample on leftmost column);
- 2. Appropriate column headings (heading for 'sample/extract/plant extract');
- 3. Observations (correct phrasing of observation with good details); R! ticks/crosses
- 4. Complete set of results + correct trend;

[4]

(e) Complete Table 1.3 to match the samples, S1, S2, S3 and S4, with each plant extract.

Table 1.3

source of plant extract	sample
a root in winter	S2
a root in spring	S4
phloem sap from the stem in summer	S3
phloem sap from the stem in winter	S1

[1]

[Total: 12]

2 Yeast contains an enzyme that will break down hydrogen peroxide into oxygen and water. The loss of mass resulting from the release of oxygen can be measured.

You will investigate the effects of the concentration of hydrogen peroxide on the rate of enzymatic activity.

You are provided with:

- Yeast cell suspension, in a container labelled Y
- 1.0% hydrogen peroxide solution, in a container labelled P
- Distilled water, in a container labelled W

Hydrogen peroxide solution P is an irritant and is harmful. Suitable eye protection should be worn. If P comes into contact with your skin, wash them off immediately under cold water.

Each time that you take a sample of yeast cell suspension Y, you should make sure that it is mixed thoroughly by stirring it with a glass rod. You should also collect Y from below the surface, so as to minimise the volume of froth collected.

(a) (i) You will carry out proportional dilutions of the 1.0% hydrogen peroxide solution P to obtain a range of concentrations in which the concentration of hydrogen peroxide is reduced by 0.2% between each successive dilution.

You will prepare 10.0 cm³ of each concentration, using **P** and **W**.

Using a table in the space below, show how you will prepare the different concentrations of hydrogen peroxide solution. One of the concentrations should include 1.0% hydrogen peroxide.

concentration of hydrogen peroxide / %	volume of P / cm ³	volume of W / cm ³
1.0	10.0	0.0
0.8	8.0	2.0
0.6	6.0	4.0
0.4	4.0	6.0
0.2	2.0	8.0

- 1. Appropriate table layout (final H₂O₂ concentration on leftmost column) + correct column heading;
- 2. 5 correct concentrations of H₂O₂ given;
- 3. Correct volumes of P and W, both to 1 dp;

[3]

7

Read steps 1–6 before starting the investigation. Proceed as follows.

- 1 Prepare the concentrations of hydrogen peroxide solution according to the table in **(a)(i)**, in the vials provided. Label the vials where appropriate.
- 2 Use the electronic mass balance to weigh the mass of the reaction mixture in subsequent steps.
- **3** Put 5.0 cm³ of yeast suspension **Y** into the vial containing 1.0% hydrogen peroxide solution. Start timing immediately and proceed to step **4** without delay.
- 4 Weigh the mass of the reaction mixture, including the weight of the vial, using the set-up shown in Fig. 2.1.



Fig. 2.1

- 5 After five minutes, weigh the mass of the reaction mixture in the vial again.
- 6 Repeat steps **3 5** with each of the other concentrations of hydrogen peroxide that you prepared in step **1**. You do **not** need to wait to complete the measurement for one concentration of hydrogen peroxide before starting on the next concentration.
 - (ii) Record your results in Table 2.1. Calculate the change in mass and the percentage change in mass. No workings are required.

		Table 2.1		
concentration	initial mass of	final mass of	change in	percentage
of hydrogen	reaction	reaction	mass / g	change in
peroxide / %	mixture / g	mixture / g		mass
1.0	20.27	20.20	- 0.07	- 0.35
0.8	20.65	20.60	- 0.05	- 0.24
0.6	20.40	20.36	- 0.04	- 0.20
0.4	20.61	20.59	- 0.02	- 0.097
0.2	20.34	20.33	- 0.01	- 0.049

- 1. All initial masses and final masses recorded to 2dp;
- 2. All change in masses correctly calculated and recorded to 2dp;
- 3. All percentage change in mass correctly calculated and recorded to 1 or 2sf;
- 4. Correct trend (highest percentage change for 1.0% or 0.8% + lowest percentage change for 0.2% or 0.4%);

[4]

may 2m

(iii) Suggest a suitable control for this experiment to show that it is an enzyme that catalyses the break down of hydrogen peroxide.

Adding 5.0 cm³ of <u>boiled and cooled</u> Y / yeast suspension to hydrogen peroxide instead; R! distilled water

.....[1]

- (iv) Other than the lack of a suitable control, describe **two** modifications to this method that would increase the confidence in your results.
 - 1. Maintain the temperature at 37°C (A! 30 to 40°C) using a <u>thermostatically</u> <u>controlled water bath</u> / using a mix of cold and hot water and measuring the temperature using a <u>thermometer</u>;
 - 2. Maintaining the pH of the reaction mixture by adding 1 cm³ of <u>buffer</u> solution;
 - 3. Repeating the experiment at least two more times and obtaining a <u>mean /</u> <u>average</u> OR to allow identification of anomalies;
 - 4. Repeating the experiment using higher concentrations of hydrogen peroxide (to obtain a more significant change in the mass of the reaction mixture);
 - **R! Use of smaller concentration intervals**
 - R! Use of vol-related apparatus with finer precision

 [2]

(b) A student wanted to determine the Michaelis-Menten constant (K_m) for the enzyme-catalysed break down of hydrogen peroxide by the enzyme found in yeast, enzyme Y.

 K_m is the substrate concentration at which the reaction rate is 50% of the maximum rate of reaction (V_{max}).

 K_m gives an indication of the affinity an enzyme has for its substrate.

The student measured the initial rate of reaction at different concentrations of hydrogen peroxide. The results are shown in Fig. 2.2.



(i) Use the graph in Fig. 2.2 to estimate the Michaelis-Menten constant (K_m).

Show your working on the graph and in the space below.

- 1. reads off V_{max} correctly from the graph at <u>4.3 au;</u>
- 2. calculates $\frac{1}{2}$ V_{max} correctly at <u>2.15 au</u>;
- 3. reads off the value for K_m (in %) correctly from the graph at <u>0.9%</u>; A! 2d.p

(ii) The K_m value for another enzyme, **Z**, is 2.5%.

State which enzyme, **Z** or **Y**, has a **lower affinity** for its substrate.

Give a reason for your answer.

enzyme Z;

reason K_m value for enzyme Z is higher / indicating a greater concentration of substrate is needed to reach a reaction rate that is 50% of V_{max} ;

......[2]

(c) Living, respiring yeast cells also contain enzymes which reduce methylene blue, a blue dye commonly used for the staining of biological samples, turning it colourless. In non-respiring cells, the reduction of methylene blue does not occur.

It is thought that respiration in yeast cells is inhibited by a high concentration of sodium chloride solution in the immediate environment.

The half maximal inhibitory concentration (IC_{50}) is a measure of the concentration of a particular inhibitory substance that is needed to inhibit a given biological process by 50%.

Plan an investigation, based on observing the colour of yeast cells mixed with methylene blue, to measure the IC_{50} value of sodium chloride solution on the respiration of yeast.

You have been provided with the following which you must use:

- prepared sodium chloride solutions with concentrations ranging from 0.1% to 2.0%
- yeast cell suspension, containing glucose as the respiratory substrate
- methylene blue solution
- microscope with an eyepiece graticule.

You may select from the following apparatus and plan to use appropriate additional apparatus:

- normal laboratory glassware, e.g. test-tubes, boiling tubes, beakers, measuring cylinders, graduated pipettes, glass rods, etc.
- syringes and Pasteur pipettes
- timer, e.g. stopwatch
- microscope slides and cover slips.

Your plan should:

- have a clear and helpful structure such that the method you use is able to be repeated by anyone using it
- identify the dependent variable and the independent variable
- identify the variables you will need to control
- use the correct technical and scientific terms
- indicate how the results will be recorded and analysed.

You can consider all steps in the procedure to be low risk and there is therefore no need to include reference to any safety measures in your plan.

.....[8]

For Examiner Use

1	Independent	State what the independent variable is and use at least five different
	variable	variables with <u>regular</u> intervals.
		The independent variable is the concentration of sodium chloride. Use
_		sodium chloride solutions of 0.4, 0.8, 1.2, 1.6 and 2.0%.
2	Dependent variable	State what is the dependent variable.
		The dependent variable is the percentage of veget calls that turn
		The <u>dependent valiable</u> is the percentage of yeast cells that turn colourloss, which is used to measure the rate of respiration
2	Donondont variable	Colouriess, which is used to measure the fate of respiration.
Т	Dependent variable	Specifies method of measuring / monitoring DV
т 1		3 Counting number of yeast cells that are colourless over the total
+		number of yeast cells
5		4. Within the field of view of the microscope using an appropriate
-		objective (need to specify how the area to be observed can be fixed)
		5. Calculate percentage of cells that turn colourless
6	Controlled	Identifies at least TWO variables to control
	variables	
		1. Volume of yeast suspension and volume of sodium chloride solution
		(to mix in test-tube/beaker)
		2. Volume of pH buffer
		3. Volume of mixture of yeast+sodium chloride and volume of
		methylene blue (to put on glass slide)
		4. Temperature of incubation of yeast+sodium chloride mixture
		5. Length of time of incubation of yeast+sodium chloride mixture
7	Controllad	6. Length of time left for staining by methylene blue
1	variables	Describes now two identified variables are controlled. (Must specify appropriate value)
+ 8	variables	appropriate value)
U		1 Using a syringe/measuring cylinder, add 5 cm ³ of yeast suspension
		to 5 cm ³ of sodium chloride solution in a test-tube / beaker:
		2. Using a syringe/measuring cylinder, add 1cm ³ of pH buffer to the
		mixture of yeast and sodium chloride solution;
		3. Using a <u>Pasteur pipette/dropper</u> , add one drop of the mixture of yeast
		and sodium chloride to one drop of methylene blue on a glass slide;
		4. Maintain the temperature at <mark>37°C</mark> (A! 30 to 40°C) using a
		thermostatically controlled water bath / using a mix of cold and hot
		water and measuring the temperature using a <u>thermometer</u> ;
		5. Start the stopwatch and let the reaction mixture incubate for 5
		minutes;
•	Decording of	6. Start the <u>stopwatch</u> and let the methylene blue stain for <u>5 minutes</u> ;
9	recording of	shows now results are to be presented in the form of a table with independent variables in appropriate columns/ rows and
	results	annronriate units for headings
		R! "Table" without lines and/ or borders
		concentration of Number of Total number of Percentage of
		sodium chloride / colourless yeast yeast cells yeast cells that
		% cells are colourless

For Examiner Use

10	Recording of results	Shows graph with appropriate axes labels and units, line/ curve with expected trend drawn
		percentage of yeast cells that are colourless
11	Calculation of IC ₅₀	Estimating the IC ₅₀ value by reading off the graph drawn where the percentage of yeast cells turning colourless is 50%.
12	Control	Describe a suitable control
		Conduct the (negative) <u>control</u> experiment using specified volume of <u>distilled water</u> instead of sodium chloride. (This shows that the inhibition of yeast cell respiration is due to sodium chloride.)
13	Reliability	Reference to repeating at least two more time with different experimental subjects
		Repeat the entire experiment at least twice using fresh yeast suspension and sodium chloride solutions, and calculate the <u>mean/average</u> of the results OR to allow identification of anomalies
14	Accuracy	Reference to relevant method to increase the accuracy of results
		Repeating experiment with smaller intervals of sodium chloride concentrations so as to obtain a better estimate of the IC_{50} .
15	Validity (other considerations)	Other considerations to improve validity of experiment
		 Mixing the solutions with a glass rod on the glass slide / AVP Before withdrawing yeast suspension, make sure to mix thoroughly by stirring it with a glass rod / collect yeast suspension from below the surface, so as to minimise the volume of froth collected. Lower coverslip gently at an angle / tap after placing the coverslip to minimize trapping of air bubbles under the coverslip

[Total: 23]

L1 is a slide of a stained transverse section through a leaf of a land plant that is affected by a fungal infection. The fungal infection affects the upper leaf surface of this plant.

(a) (i) Use your microscope to observe the different tissues in the region of slide L1 shown by the darkly shaded area in Fig. 3.1. The observed area should include at least one vascular bundle and be affected by fungal infection.



Fig. 3.1

Draw a large plan diagram of the part of the leaf shown by the shaded area in Fig. 3.1.

A plan diagram shows the arrangement of different tissues, including their correct shapes and proportions. No cells should be drawn.

Your drawing should also include any fungal tissue observed. Use **one** ruled label line and label to identify the fungal tissue.



1.	Accuracy	Show all the structures that can be seen in the defined part (2 epidermal layers + 2 mesophyll layers + vascular bundle)			
2.	Accuracy	Correct shape and proportions of tissue types (spongy mesophyll layer slightly thicker than palisade mesophyll layer) <i>R! individual cells/ individual xylem vessel drawn</i>			
3.	Clarity	Use of sharp pencil to draw clear, single, continuous, neat lines			
4.	Scale	Use at least 2/3 of space provided			
5.	Label	 a) Fungal tissue of correct proportion drawn at the upper epidermal layer b) Clear, straight line with label for 'fungal tissue' 			



[5]

(ii) Observe one vascular bundle of the section on L1.

Select one large xylem vessel element and three adjacent smaller cells.

Each smaller cell must touch the large xylem vessel element and at least one of the other smaller cells.

Make a large, labelled drawing of these four cells.



1.	Accuracy	a) 4 irregularly shaped cells (1 large + 3 smaller cells)					
		b) Each smaller cell touching the large cell and one					
		other cell					
2.	Accuracy	Shows thick cell wall in at least the large cell					
		R! chloroplasts drawn in smaller cells					
3.	Clarity	Use of sharp pencil to draw clear, single, continuous, neat					
		lines;					
4.	Scale	Use at least 2/3 of space provided					
5.	Label	a) Clear straight label lines					
		b) Labels for: (thickened) cell wall					
		R! cytoplasm label for large xylem vessel					

[Turn over



[5]

(iii) Fig. 3.2 is a photomicrograph of a stained transverse section of part of a leaf from a different species of plant.

Observe the photomicrograph in Fig. 3.2 and the section on **L1** to identify differences between them.

Fig. 3.2 has been annotated to describe **one** of these differences. A label line has been used to indicate the feature that is different.

Complete Fig. 3.2 by:

 identifying and annotating two more differences between the micrograph in Fig. 3.2 and the section on L1

Fig. 3.2 Upper epidermis is continuous and not disrupted

• using a label line to identify the feature that is different.

L1 Upper epidermis is disrupted by fungal tissue

Fig 3.2

Any of the two comparisons between Fig. 3.2 vs L1;;

- 1. Mesophyll layers: No clear distinction of palisade and spongy mesophyll layers vs clear upper layer of palisade mesophyll and lower layer of spongy mesophyll;
- 2. Arrangement (or number) of vascular bundles: Vascular bundle arranged in two/three rows vs vascular bundle arranged in a single row;
- 3. Size of vascular bundles: Vascular bundle is small relative to the leaf section/does not span the leaf section vs vascular bundle is large relative to the leaf section/may span the leaf section;
- 4. Arrangement (or number) of air spaces: Two/three rows of air spaces vs single row of air spaces;
- 5. Size of air spaces: Air spaces are large relative to the leaf section/may span the leaf section vs air spaces are small relative to the leaf section/does not span the leaf section;
- 6. Label line correctly drawn to feature being compared; *compulsory MP

[3]

[Turn over

(b) A scientist investigated changes in the mean width of stomata in the leaves of a plant growing in hot, dry conditions. The scientist measured the widths of stomata at different times of day, from 02:00 hours to 22:00 hours. Fig. 3.3 shows where the scientist measured the width of each stoma.



Fig. 3.3

The scientist calculated the mean width of stomata for each time of day. The results are shown in Table 3.1.

time of day	mean width of stomata		
/ hours	/ arbitrary units (au)		
02:00	86		
04:00	36		
07:00	4		
15:00	2		
22:00	95		

т	a	b	le	3.	1
	~	~		•	



- 2. Axes scaled appropriately so that graph takes up at least 50% of the grid on both axes and divisions are equidistant; A! intervals of 4h or 5h on x-axis R! zero time represented as '0' instead of '00:00' R! 25:00 if using 5h on x-axis
- 3. All points correctly plotted;
- 4. Dot-to-dot plot or line of best fit, with no extrapolation beyond extreme measured data;

[4]

For

Use

(ii) Fig. 3.5 shows a stage micrometer scale that is being used to calibrate an eyepiece graticule.

The length of one division on this stage micrometer is **0.02 mm**.



Using the eyepiece graticule shown in Fig 3.5, the width of a guard cell measures **3 eyepiece graticule divisions**.

Calculate the actual width, in micrometres (µm), of this guard cell.

Show all the steps in your calculation, including the appropriate units.

- Calibrating eyepiece graticule using stage micrometer:
 5 stage micrometer divisions = 50 eyepiece graticule divisions = 0.02 x 5 mm
 A! ratio of 1 stage micrometer division : 10 eyepiece graticule divisions
- 2. Calculating length of one eyepiece graticule division: 1 eyepiece graticule division = $0.02 \times 5 \div 50 \text{ mm} = 0.002 \text{ mm} = 2 \mu \text{m}$
- 3. Calculating actual width of guard cell + conversion into μm (can be done in earlier steps):
 3 eyepiece graticule divisions = 3 x 2 μm = 6 μm (1sf) A! 2sf

R! if no description given for working step

actual width of guard cell = µm

[3]

[Total: 20]